## **GROUP B CATIONS**

Mn <sup>2+</sup>	Fe <sup>3+</sup>	Bi <sup>3+</sup>	Cr <sup>3+</sup>
	. •	<b>_</b> .	<u> </u>

Group B cations are selectively precipitated out using ammonia buffer at pH 9-10 to form hydroxides. Ammonia provides the right amount of hydroxides ions in the solution to selectively precipitate group B cations, and forms stable complexes with the transition metal ions to make sure that these remain in soluble form in the buffered solution:

$$\begin{array}{lll} Bi^{3+}_{(aq)} \ + \ NH_4OH_{(aq)} \ \rightleftharpoons \ Bi(OH)_{3(s)} + NH_4^+_{(aq)} \\ \\ Fe^{3+}_{(aq)} + NH_4OH_{(aq)} \ \rightleftharpoons \ Fe(OH)_{3(s)} + NH_4^+_{(aq)} \\ \\ Mn^{2+}_{(aq)} + NH_4OH_{(aq)} \ \rightleftharpoons \ Mn(OH)_{2(s)} + NH_4^+_{(aq)} \\ \\ Cr^{3+}_{(aq)} + NH_4OH_{(aq)} \ \rightleftharpoons \ Cr(OH)_{3(s)} + NH_4^+_{(aq)} \end{array}$$

Mn <sup>2+</sup>	<ul> <li>Confirmatory reagent: 3% H<sub>2</sub>O<sub>2</sub>/NaBiO<sub>3(s)</sub></li> <li>Reduction of Mn<sup>4+</sup> to Mn<sup>2+</sup> by reacting with H<sub>2</sub>O<sub>2</sub> first, allows the reaction with BiO<sub>3</sub><sup>-</sup> to occur</li> <li>Production of purple supernatant</li> </ul>	
Fe <sup>3+</sup>	Confirmatory reagent: 0.1 M KSCN	
	<ul> <li>Formation of blood red thiocyanatoiron (III) complex</li> </ul>	
Bi <sup>3+</sup>	Confirmatory reagent: 6 M NaOH/SnCl <sub>2(s)</sub>	
	<ul> <li>Bi<sup>3+</sup> is reduced to Bi with Sn<sup>2+</sup> as the reducing agent</li> </ul>	
	<ul> <li>This reduction reaction occurring in a basic solution allows</li> </ul>	
	the immediate appearance of a black precipitate	
Cr <sup>3+</sup>	Addition of H <sub>2</sub> O <sub>2</sub> + HCl	
	- Reduces Cr₂O²- to CrO₅	
	- A flash of <mark>blue</mark> color in the solution, which decomposes	
	quickly, is observed	