

GENERAL FLOWCHART FOR GROUP A CATIONS

**ORIGINAL SOLUTION:
SOURCE OF CATIONS**

Test your sample for precipitation and centrifuge

RESIDUE

CENTRIFUGATE

Wash sample (solids) with 10 drops of **COLD** distilled water. Add 2 drops of room temperature distilled water.

Heat in boiling water bath for 3 minutes with constant stirring. Centrifuge quickly and decant supernatant liquid while sample is **HOT**.

**CA1: CATIONS
FOR OTHER
GROUPS**

RESIDUE

CENTRIFUGATE

Wash it with 10 drops of **HOT** distilled water. Add 2 drops of 6M NH_3 . Stir and Centrifuge

Add 3 drops of 1 M KI

The sample should produce **YELLOW** precipitates (solids) to indicate presence of PbI_2 .

RESIDUE

CENTRIFUGATE

Your residue appears to be **BLACK GRAY** or **BLACK**

Add 6 M HNO_3 dropwise with stirring until acidic
The **blue** litmus paper should turn **RED** which means it is now acidic.

Your sample is Pb^{2+}

Your sample is Hg_2^{2+}

Your sample produced **WHITE** precipitates/solids

Your sample is Ag^+